



# E-xtreme

Make Essence to Extreme  
Smaller | Lighter | Quieter  
than ever


# Smaller



• 360° contra-angle rotation



• Miniature contra-angle for better visibility & easier posterior access




Miniature contra-angle  
Super Mini head & thin neck



Normal head & neck


• Slim and compact design for maximum comfort and mobility



# Lighter



Less than 100 grams,  
reducing hand fatigue during treatment,  
providing the ultimate E-xtreme experience.



# Quieter

Eighteeth innovative motor,  
cutting down the noise from inside,  
allowing you to experience the tranquility.



Reciprocation function:  
five optional angles



Built-in file system



Model	: E-xtreme
Weight(grams)	: 99 g
Torque range	: 0.5 N.cm-4.0 N.cm
Power supply	: 100 V-240 V
Frequency	: 50 Hz / 60 Hz
Output	: DC 5 V/ 1 A