



WIRELESS Z

Wireless Headlight

Wireless Light, Clarity Get



29.5_g
Super Light



600_{mAh}
Powerful Battery

Touchscreen Brightness Control



124_{mm}
Light Spot at
450mmWD



Up to
35000_{lux}


CRI≥90%

Wireless



Wireless Headlight

WIRELESS LIGHT, CLARITY GET



Wireless

Cable-free stable lithium batteries ensure a smooth work experience.



Battery Status Indicator

Built-in battery status indicator flashes when the battery is low.



Brightness Memory

Automatically adjust to brightness when last powered off to reduce adjustment time and improve work efficiency.



Touch Control Design and Stepless Dimming

Sensitive touch technology can be used with gloves.



Super Light

Weigh just 27g (battery, headlight, filter included) and sit on the frame's upper edge for reducing nose bridge pressure.



Powerful battery

Ensure Uninterrupted Surgical Illumination



Long-lasting


Sustain maximum brightness for up to 2.5 hours per battery.
Package comes with two rechargeable batteries.

Fast

Amazing charging speed can greatly save your time.
(Single charging for 1.4 hours and double batteries for 2.6 hours.)

Convenient


Built-in battery status indicator. Magnetic positioning for swift and easy charging.



124mm

Light Spot at 450mmWD

Providing clear and uniform illumination that effectively lights the entire arch with new optical design.



Up to 35,000 LUX

Even edge-to-edge and same light intensity spot to ensure clarity.



CRI $\geq$ 90%

Restoring true colors with industry-leading imported medical lighting LEDs.



Compatible with Eighteeth® Brilliance Series Loupes.

(Optional: customized connections with specified other brand loupes.)



Magnetic connection
for quick attachment.

SPECIFICATIONS

Total Weight: 27g

Battery Capacity: 550mAh, 3.7V

Charging Time: Single charging 1.4h; Dual charging -2.6h


Battery Type: Li-ion

Adapter Input: 100-240 V 50/60 Hz, 0.2 A

Adapter Output: 5V1A

Lamp weight(including adaptor): 12g

Battery Module: 15g



Eighteeth